

# Advanced abdomen for injectables practice

PRODUCT NUMBER: 1436



## DESCRIPTION:

The Advanced Abdomen for Injectable Practice is an invaluable tool for demonstrating how to correctly insert and rotate the insulin pump infusion set, avoiding the 5 cm zone around the umbilicus, as well as for administering many types of self-injections (including pen injections).

## CHARACTERISTICS:

### GENERAL:

- Made of high quality material
- Realistic skin with what feels like real stomach tissue
- Ability to pinch the skin
- Soft material can be repeatedly punctured with no visible holes. Needle resistant to needle through
- Adjustable straps allow the stomach to be attached to a mannequin or human patient for simulation practice
- Fluid injection is not recommended
- Dimensions: 25.5 x 16 x 5cm
- Manufactured under an ISO 9001 certified quality management system

## INCLUDE:

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

- 1 Advanced Abdomen Model for injectables practice
- 1 Adjustable strap



**MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

**CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>