



PERICARDIOCENTESIS SIMULATOR

PRODUCT NUMBER: 2296



DESCRIPTION:

The **Pericardiocentesis Simulator** is designed specifically to teach the skills needed to perform this difficult procedure correctly, as well as ongoing chest tube maintenance and the **management of pre-hospital chest trauma**.

SKILLS:

- Chest tube surgical placement and maintenance
- Pericardiocentesis procedure skills training
- Complete pre-hospital trauma management
- Needle decompression of pressurized pneumothorax

CHARACTERISTICS:

GENERAL:

- The simulator has a pressurized tension pneumothorax site and a site for the surgical placement of a functional chest tube.
- Fluid color, volume, and viscosity can be controlled by the instructor.
- The manikin may be used with any commercially available closed water seal drainage unit.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Pericardiocentesis may be performed in the left subxiphoid space and in the left fifth intercoastal space.
- Accurate placement of the needle will allow for the withdrawal of fluid from the simulated pericardium.

INCLUDE:

- Pericardiocentesis simulator
- 1 Simulated pericardium
- 1 Fluid reservoir bag
- 1 Foot pump
- Surgical skin pads
- Subcutaneous surgical pads
- Nurse training pads
- Pneumothorax pads
- Simulated blood concentrate
- 1 Hard carry case

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>