

# MULTIPURPOSE VENOUS TRAINING ARM

PRODUCT NUMBER: 1165



## DESCRIPTION:

The Multipurpose Venous Training Arm provides the ability to train, teach and practice intramuscular and subcutaneous injections with all the features of an intravenous (IV) arm. Train and teach students on medication administration simulation, blood collection practice, and aseptic techniques in a safe, sanitary learning environment.

\* Vein tubes contain latex, which may cause allergic reactions. Users allergic or sensitive to latex should avoid contact. Discontinue use of this product and seek medical attention if an allergic reaction occurs.

## SKILLS:

- IV training exercises
- Blood collection
- Infusion exercises
- Intravenous injections
- Intramuscular injections
- Subcutaneous injections
- Venipuncture
- Toruniquet technique
- Venous therapy
- Detection of collapsed or bulging vessels

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Professional-to-patient communication

## CHARACTERISTICS:

### ANATOMY:

- Subtle multi-venous network
- Realistic “pop” as needle enters vein
- Resealing veins and outer skin

### VERSATILITY:

- Subcutaneous injection areas on the volar side of the forearm and the lateral side of the upper arm
- Intramuscular injection site in deltoid area
- TB test site
- Squeeze bulb to increase or decrease venous pressure. Veins stand out or collapse.
- Administration of medication by intravenous bolus.
- Simulation of infusion technique
- Blood collection exercises with simulated blood
- Simulation of clenched fist and tourniquet position

### USER-FRIENDLY:

- Easily replaceable skin and veins
- Compact, thereby improving portability of the unit.
- Easy to clean, service and maintain

## INCLUDE:

- 1 Multipurpose IV training arm
- 1 Smoked Lucite base with stand for blood bag
- 1 Blood dispensing bag
- 1 Pressure bulb
- 1 Spare arm skin
- 1 Carrying bag
- 1 Instruction manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>