



# NEWBORN INJECTION TRAINING LEG

PRODUCT NUMBER: 1171



## DESCRIPTION:

The **Susie® and Simon® intraosseous leg** simulates that of a newborn.

A **replaceable tibial bone lies under a smooth outer skin**. It includes anatomic landmarks for teaching intraosseous access and infusion. A femoral vein/artery pair and an intramuscular injection site are included.

## SKILLS:

- Intraosseous access
- Infusion of medicaments/fluid exercises
- Blood draw exercises
- Intramuscular injection
- IV training exercises

## CHARACTERISTICS:

- Newborn-sized right leg
- Mounted on smoked Lucite® base

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Blood bag and stand provide continual blood supply
- Intraosseous vascular access
- Realistic “pop” as needle penetrates cortex and enters bone marrow
- Blood return upon successful entry
- Veins withstand multiple needle sticks
- Easily replaceable skins and veins
- Sixteen interchangeable, realistic tibial bone inserts

## INCLUDE:

- 1 Newborn IV training leg
- 1 Femoral vein
- 1 Stand
- 16 Interchangeable tibial bones
- 1 Blood bag
- 1 Squeeze bulb
- 1 Soft carrying bag
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>