

# HAL EMERGENCY AIRWAY & CPR - OMNI2

PRODUCT NUMBER: 1223



## DESCRIPTION:

HAL® S315.100.250 is an easy-to-use portable airway management and CPR skills trainer designed for teaching resuscitation and emergency airway management skills to students and professionals in EMS, Medicine, and Nursing. HAL® features a high-fidelity airway that simulates laryngospasms, true-to-life chest recoil, and real-time CPR quality feedback to help you teach effective basic and advanced airway management as well as cardiopulmonary resuscitation.

## SKILLS:

- Airway management
- Nasal and oral intubation
- Resuscitation training (CPR)
- BVM ventilation

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Emergency airway management (tracheotomy or cricothyrotomy)
- Detection of laryngospasm
- Debriefing tool

## CHARACTERISTICS:

### GENERAL:

- Adult-sized, male upper torso
- Anatomically correct landmarks for proper CPR hand placement

### AIRWAY MANAGEMENT:

- Realistic airway with tongue, vocal cords, and trachea
- Realistic head tilt/chin lift and jaw thrust
- Nasal and oral Intubation with standard adjuncts: ETT, LMA, King LT®
- CPR hand position landmarks
- Visible chest rise with bag-valve-mask ventilation
- Esophageal intubation and/or excessive ventilation produce gastric distension
- Laryngospasms
- Supports tracheostomies and cricothyrotomies
- Right mainstem intubation presents unilateral chest rise

### CPR:

- Compliant with 2015 adult CPR guidelines
- Realistic chest cavity resistance and recoil
- Real-time feedback of CPR quality: Compression depth and rate, Ventilation Rate, Excessive Ventilation, No-Flow Time and CPR Cycles

## INCLUDE:

- 1 HAL Skills trainer
- 1 OMNI2 Wireless Control Interface (Tablet+License)
- 1 Carrying case
- 1 Power supply adapter
- 1 Instructions manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>