

BREATHLESS PREMIE SUZIE SKILLS TRAINER

PRODUCT NUMBER: 2059



DESCRIPTION:

Suzie is a freestanding model of a premature born baby girl with breathing problems. She has an open mouth which shows a blockage.

SKILLS:

Tracheotomy skills training with premature patients

CHARACTERISTICS:

- The easy replaceable neck skin can be penetrated and under the muscles tissue the trachea and larynx is visible.
- There are 2 bloodlines on either side of the trachea that can be connected to the syringe for lubrication (diluted shampoo or liquid soap) or blood if that is needed.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Suzie is made with soft silicon and hand painted with rooted hair and eyelashes.
- This model is easy to use many times over.
- The neckskin can be repaired or replaced.

INCLUDE:

- 1 Baby Suzie
- 5 Neck skins
- 3 Trachea inserts
- 1 Syringes

