

# PEDIATRIC INJECTION TRAINING - 5 YEAR

PRODUCT NUMBER: 1167



## DESCRIPTION:

The Mike® and Michelle® injection training arm and hand simulates that of a 5 year old, with lifelike venous structure. The dorsum of the hand has two veins for additional intravenous set-up and training. The flexor and lateral sides of the arm contain sites for subcutaneous injection practice. The deltoid contains an intramuscular injection site.

## SKILLS:

- IV training exercises
- Venipuncture training
- Blood draw exercises
- Infusion exercises
- Intramuscular injections
- Subcutaneous injections
- Venous therapy
- Tourniquet technique

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Detection of collapses or bulging vessels
- Professional-to-patient communication

## CHARACTERISTICS:

### FEATURES:

- VERSATILITY
  - Intravenous injections, infusion, and blood collection training
  - Subcutaneous injections in both flexor and lateral sides of forearm
  - Intramuscular injections in deltoid
- ANATOMY
  - Lifelike venous structure
  - Dorsum of hand has two veins for additional intravenous set-up
  - Squeeze bulb simulating either collapsed or distended veins
- USER-FRIENDLY
  - Easily replaceable skin and veins that withstand multiple needle sticks
  - Compact, thereby improving portability of the unit.

## INCLUDE:

- 1 Pediatric IV training arm
- 11 Squeeze bulb with 2.5ft line
- 1 Smoked Lucite Base with stand for blood bag
- 1 Spare arm skin
- 1 Dispensing bag
- 1 Synthetic blood concentrate
- 1 Carrying bag
- 1 Instruction manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>