

PEDI BLUE Newborn CPR patient simulator (SmartSkin) - OMNI1

PRODUCT NUMBER: 1239



DESCRIPTION:

The PEDI® Blue Neonatal manikin is a spectacular neonate newborn simulator which changes cyanosis color based upon an initial pre-selected condition and measures the effectiveness of CPR, airway ventilation and chest compression. In addition, the simulator has all the conventional features found in airway management trainers.

Optional accessories include an intraosseous leg and an injection training arm.

SKILLS:

- Pediatric airway management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Resuscitation skills (CPR)
- BVM ventilation
- Nasal and oral intubation
- Umbilical catheterization
- Cyanosis detection and management
- Intraosseous infusion (with optional leg)
- Injection and intravenous exercises (with optional arm)

CHARACTERISTICS:

FEATURES:

- Prominent venous network
- Full size articulating neonate
- Realistic airway with tongue, vocal cords, trachea, and esophagus for airway management exercises
- Internal organs for unparalleled CPR performance
- BVM with Realistic Chest Rise Bilateral lung expansion with realistic chest rise during BVM.
- Oral or nasal intubation plus suctioning
- Crico prominence permits the Sellick maneuver
- Bilateral lung expansion with realistic chest rise
- Practice intubation using a Miller 1 blade and uncuffed 2.5 ET tube or LMA
- View peripheral and central cyanosis as well as healthy skin tone
- Practice umbilical catheterization and simulate pulsatile element with use of hand-held squeeze bulb
- Simulate neonatal CPR with either two-thumb "encircling" technique or two finger alternate compression method
- Practice delicate IO access with optional intraosseous leg (sold separately)
- Practice injection and intravenous techniques with optional training arm (sold separately)

INCLUDE:

- 1 PEDI Blue neonate simulator
- 1 OMNI1 controller
- 1 Power supply adapter
- 1 Carrying case
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>