


Nursing care simulator - 5 year

PRODUCT NUMBER: 1271



DESCRIPTION:

The Five-Year-Old Multipurpose Patient Simulator is designed to help participants practice and develop pediatric nursing care and resuscitation skills, including CPR, intubation, trach care, IV access, and more.

SKILLS:

- Resuscitation skills (CPR)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- BVM ventilation
- Nasal and oral intubation
- Nasogastric intubation
- Gastric suctioning exercises
- Male and female urinary catheterization
- Tracheostomy care
- Intraosseous exercises
- Intravenous exercises
- Intramuscular, subcutaneous and intradermal sites for injection
- Enema administration
- Estoma care and management
- Placement of rectar suppositories
- Pulse palpation
- Ear canal examination
- Administration of orbital medicines
- Removal of foreign bodies
- Eye irrigation
- Bandaging

CHARACTERISTICS:

FEATURES:

- Full-body with soft, lifelike face skin and molded hair
- Eyes can be closed or opened manually. It has a dilated left pupil for advanced exercises.
- Articulating head, jaw, arm, and legs
- Interchangeable genitalia for female and male catheterization
- Lavage/gavage
- G-tube placement
- Enema administration
- Placement of rectal suppositories
- Simulated ear canal
- Practice gastric suctioning exercises

AIRWAY MANAGEMENT:

- Endotracheal intubation
- Head tilt/chin lift
- Nasogastric intubation
- Insert a tracheostomy tube with an attached positive pressure device to observe chest rise.
- Tongue, teeth, and realistic airway

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

CPR:

- CPR anatomical landmarks
- Accommodates assisted ventilation
- BVM with realistic chest rise

VENOUS ACCESS:

- Intraosseous infusion and injection system with realistic tibia bones.
- IV training arm and hand for intravenous subcutaneous and intradermal exercises
- Femoral venous access
- Intramuscular Injection sites on the deltoids and quadriceps

CIRCULATION:

- Manual palpable pulses: carotid, femoral, brachial
- Stomas for ileostomy, colostomy, and suprapubic exercises

OPHTHALMIC EXERCISES:

- Administration of orbital medicines into the conjunctival sac
- Removal of foreign bodies
- Eye irrigation

INCLUDE:

- 1 Pediatric patient simulator manikin
- 1 Interchangeable male genitalia
- 3 Stomas (installed)
- 2 Filling syring
- 5 IO tibia bone
- 2 Filling syring
- 1 Carrying bag
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>