



ADVANCED CPR DOG

PRODUCT NUMBER: 1953



DESCRIPTION:

The advanced sanitary CPR dog provides practice in assisted breathing, students may practice the appropriate rate, pressure, and position of cardiac massage; and coordinate the respiratory cardiovascular functions.

SKILLS:

- CPR training
- Ventilation training
- IV catheter insertion

CHARACTERISTICS:

- The advanced sanitary CPR dog provides each student with their own nosepiece with a one-way valve. The nosepieces and lower airways are disposable.
- In addition to practice in assisted breathing, students may practice the appropriate rate, pressure, and position of cardiac massage; and coordinate the respiratory cardiovascular functions.
- The Advanced CPR Dog is connected to a small box that displays a light when the correct ventilation volume is applied, another light

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

indicates correct position for cardiac massage, a third light indicates correct compression depth, and an audible beep and light indicates excessive pressure is being applied.

- A visible chest rise appears when correct ventilations are applied.
- An IV can even be established for more advanced training.

INCLUDE:

- 1 CPR trainer dog
- 1 Soft carrying case
- 10 Nosepieces
- 3 Lower disposable airway.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>