



Male and Female Catheterization Trainers

PRODUCT NUMBER: 1331



DESCRIPTION:

Practice male and female urinary catheterization on Simulaids' newly designed anatomically correct models. Training is enhanced by a new material that has a realistic feel and flexibility.

CHARACTERISTICS:

GENERAL:

- Both models utilize a 16 FR catheter and can be used to demonstrate balloon inflation procedures.
- The models are individually mounted on a training stand and can be used on a tabletop or for patient simulation.
- A fluid reservoir allows for urine output upon successful completion of the insertion procedure.
- Instructors are able to teach sterile technique of insertion and removal of catheter, as well as practice preparing the patient for the procedure.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Includes both models listed above.
- Dimensiones: aprox. 15cm
- Weight: 3kg

INCLUDE:

- 1 Female catheterization trainer
- 1 Male catheterization trainer