


Pad for IV injection in large animals

PRODUCT NUMBER: 3361



DESCRIPTION:

This realistic looking intravenous injection pad offers the opportunity to practice blood sampling and venous catheterisation on large animals.

This blood collection pad has two integrated venae cavae, which perfectly mimic the characteristics of subcutaneous blood vessels.

SKILLS:

- The two veins simulate different levels of blood stasis and therefore allow training at different levels of difficulty, as is often the case in real life.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- The skin surface, which resembles that of a large animal, has a very realistic texture and feel.
- The skin-like surface closes after injection and allows for multiple IV injections. Reusable
- The trainer has a holder that ensures a secure hold on the table with suction cups.
- Tear resistant
- Retains its shape and texture even after repeated use.
- The holder can be used to mount the injection pad vertically or horizontally as required.
- Dimensions: 15 x 16 cm
- Thickness: 3 cm
- Veins: 0.9 cm ?

INCLUDE:

- 1 Pad for IV injection in large animals
- 1 Stand

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>