

CHESTER CHEST - Upper torso for vascular access training

PRODUCT NUMBER: 1321



DESCRIPTION:

Chester Chest™ has been an industry standard in healthcare simulation since 1986 for the teaching of central line care. Chester Chest™ features a standard arm with a dual lumen 5FR PICC exiting the basilic vein from the inner bicep area and a pre-positioned 20g IV catheter in the forearm.

SKILLS:

- Teaching of central line care.
- Practice infusion, withdrawal, securement, and dressings of multiple vascular access lines.
- IVAD placement in the chest

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Jugular and subclavian catheter placement
- Blood sampling
- Administration of medication via infusion

CHARACTERISTICS:

- Chester Chest™ makes it possible to teach and practice the infusion, withdrawal, care, securement and dressings of the following vascular access lines: IVAD placed in the chest, Hickman-like central line, PICC line, IV site and optional Triple lumen jugular and subclavian catheters.
- With this model, it is possible to simulate accessing the following IVAD placements: normal, "tipping", "wandering", or deeply placed.
- Successful access is confirmed by a blood return.
- Fluid can be infused and blood withdrawn from all the lines.
- External jugular or subclavian catheters can be added through existing openings on the model.
- Chester Chest can be used in either an upright or supine position.
- Chester Chest is a complete unit consisting of a life-sized torso with neck and chin, detachable right arm, outer tissue ap, three difficult accessing inserts, an actual implanted port, CVC, PICC, openings for attachment of a subclavian or external jugular catheter, "blood" reservoir bags, tubing, talc baby powder and user's information.

INCLUDE:

- 1 Chester Chest simulator
- 1 Chester Chest standard arm
- 1 Blood reservoir bag
- 1 Talcum powder
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>