

# ADULT X-RAY DENTAL MANIKIN

PRODUCT NUMBER: 3044



## DESCRIPTION:

Provides simulation of dental radiography adult patient preparation.

## CHARACTERISTICS:

### FEATURES:

- Includes special DENTOFORM model with radio-opaque plastic teeth, flexible finger for holding film, bite-opening instruments, and latex tongue.
- Includes chair mount
- Weight: 2,75kg

## INCLUDE:

- 1 Adult X-ray dental manikin
- 1 Chair mount

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>