

PEDI BLUE - Newborn ALS patient simulator (SmartSkin) - OMNI2

PRODUCT NUMBER: 1236



DESCRIPTION:

A spectacular newborn simulator which changes color based upon an initial pre-selected condition and measures the effectiveness of airway ventilation and chest compression. In addition, the simulator has all the conventional features found in airway management trainers.

SKILLS:

- Pediatric airway management
- Resuscitation skills (CPR)
- BVM ventilation
- Nasal and oral intubation

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- NG exercises
- Umbilical catheterization
- Intraosseous infusion
- Injection and intravenous exercises
- Palpation of pulses
- Lavage/gavage
- Debriefing exercises

CHARACTERISTICS:

FEATURES:

- Full-body, full-term infant with articulating head, arms, legs
- Realistic airway with tongue, vocal cords, trachea, and esophagus for airway management exercises
- Support oral or nasal intubation plus suctioning using standard devices
- Realistic lung expansion with bag valve mask ventilation
- Programmable peripheral and central cyanosis as well as healthy skin tone
- Intravenous (IV) access on right arm and lower left leg support infusion and bolus injections
- Right tibia intraosseous access and infusion
- Palpable pulses: right brachial, radial, femoral and umbilical (manual bulb)
- Practice umbilical catheterization and infusion
- Simulate neonatal CPR with either two-thumb "encircling" technique or two finger alternate compression method
- Supports gastrostomy tube placement

INCLUDE:

- 1 PEDI Blue neonate simulator
- 1 OMNI2 Wireless interface (tablet+license)
- 2 Replacements for tibia IO access
- Power supply adapter
- 1 Carrying case
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>