


# Pediatric Simulator PEDI (5 years)

PRODUCT NUMBER: 1273



## DESCRIPTION:

The Five-Year-Old PEDI® Patient Simulator is designed to help participants practice and develop pediatric nursing care and resuscitation skills. The S155 includes an intubatable airway, an injection training arm, an intraosseous injection leg, and appropriate arterial/venous system to ensure participants in ALS training programs receive hands-on training.

## SKILLS:

- Resuscitation skills (CPR)
- BVM ventilation
- Nasal and oral intubation
- Nasogastric intubation
- Male and female urinary catheterization

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Tracheostomy care
- Intraosseous exercises
- Intravenous exercises
- Intramuscular, subcutaneous and intradermal sites for injection
- Pulse palpation
- Bandaging

## CHARACTERISTICS:

### FEATURES:

- Full-body simulator with soft, lifelike face skin and molded hair
- Eyes can be closed or opened manually
- Head tilt/chin lift
- Endotracheal intubation
- Nasogastric intubation
- Tongue, teeth, and realistic airway
- Intraosseous infusion
- IV training arm
- Femoral vein
- Manual palpable pulses: carotid, femoral, brachial
- BVM with realistic chest rise
- CPR anatomical landmarks
- Lavage/gavage
- Intramuscular injection sites
- Simulated ear canal
- Practice gastric suctioning exercises
- Articulating head, jaw, arm, and legs

## INCLUDE:

- 1 Pediatric patient simulator manikin
- 1 Interchangeable male genitalia
- 5 IO tibia bone
- 2 Filling syring
- 1 Carrying bag
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>