

AirSim Child Bronchi X

PRODUCT NUMBER: 1104



DESCRIPTION:

TruCorp AirSim Child Bronchi X provides pediatric doctors, pediatric clinical specialists and emergency responders a way to practice pediatric airway management skills and diagnostic techniques on an anatomically correct manikin. Suitable for training in laryngoscopy, bronchoscopy, intubation and lung isolation techniques. Features the AirSim X airway with 5-year warranty.

SKILLS:

- Oral and nasal intubation
- Bag/mask ventilation
- Laryngeal Mask Airway (LMA) and supraglottic device insertion
- Combitube insertion; Double lumen endotracheal tube insertion
- Diagnostic bronchoscopy

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Lung isolation techniques using left and right endobronchial tubes and bronchial blockers
- Lung suctioning techniques

CHARACTERISTICS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- Features bronchi which provides internal anatomically correct features down to the fourth generation
- Inflatable tongue bulb to create tongue edema providing a variance during training
- User feedback through the visualization of stomach inflation and lung expansion
- Missing front tooth to illustrate realistic dentals of a 6-year old
- Modular design allowing for all parts to be changed within your organization (if required)
- Delivered ready to use in a durable carry case for easy transportation and safe storage

INCLUDE:

- 1 AirSimChildBronchi X (assembled and ready to use)
- 1 TruCorp carry case
- 1 USB pen drive containing product user manual
- 1 bottle of TruCorp lubrication (100ml)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>