



GOLDIE K9 - CANINE AUSCULTATION SIMULATOR

PRODUCT NUMBER: 2468



DESCRIPTION:

The highlight of this breath/heart sound (BHS) simulator is to allow students to perform auscultation, offering them opportunities to hear the different sounds with a stethoscope.

SKILLS:

- Canine heart sounds auscultation
- Canine lung sounds auscultation

CASES AND PATHOLOGIES:

HEART SOUNDS MODULE:

- Atrial Fib
- Mitral Regurgitation
- Mitral Valve Click
- Normal Heartbeat
- PDA
- Pulmoic Stenosis

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Respiratory Crackles and MR Murmurs
- SAS
- Mitral Regurgitation
- VPC
- VSD

LUNG SOUNDS MODULE:

- Tracheal
- Vesicular
- Bronchial-Vesicular
- Wheezes
- Monophonic Wheeze
- Pleural Friction Rub
- Stridor
- Cavernous
- Crackles
- Pulmonary Edema
- Puppy

CHARACTERISTICS:

- Instructors can choose various sound scenarios for the classroom.
- This product uses plug-in modules, with data collated from actual patients.
- Five speakers are installed directly in our Goldie canine BHS mannikin.
- The BHS unit also has two lights that illuminate during expiration. It is powered by a 9V battery and comes with a limited warranty.
- Each Goldie comes with a tutor MS, one Breath Sound module, one Heart Sound module, five installed custom speakers and a 9V battery.
- A rugged transport case with foam insert for the Breath/Heart Sounds module can be ordered as an option

INCLUDE:

- 1 Canine auscultation trainer
- 1 Breath Sound module
- 1 Heart Sound module
- 5 Installed custom speakers
- 1 Battery (9V)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>