

# CASPER CPR DOG

PRODUCT NUMBER: 2050



## DESCRIPTION:

CasPeR the CPR Dog provides training in mouth-to-snout resuscitation with individual-use, disposable airway bags.

## SKILLS:

- Proper and hand placement for compressions
- Airway management
- CPR
- Ventilation
- Femoral pulse check

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



## CHARACTERISTICS:

- These economical disposable shield/lung bags eliminate cross contamination between students and the need to disinfect the manikin after every use.
- This dog manikin can be used to teach proper hand placement for compressions, depth of compressions, airway maintenance, femoral pulse check, and adequate ventilation.
- The model has no indication function, the teacher has to judge proper exercise.
- Airways are changed quickly and without tools.

## INCLUDE:

- 1 Casper CPR dog
- 100 disposable airways

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>