

CPR Simulator with Indicator Light - Adults, Children and Infants

PRODUCT NUMBER: 2968



DESCRIPTION:

The PRESTAN Professional Collection combines two (2) professional adult mannequins, one professional child and two (2) professional child mannequins in one collection. Fully compliant with current industry guidelines, including the latest AHA Integrated Feedback Directive.

Available in medium or dark skin tones and CPR only.

CHARACTERISTICS:

GENERAL:

- All PRESTAN professional manikins are available with a revolutionary CPR rate monitor that allows instant feedback for both instructors and students.
- The visual feedback from the CPR rate monitor gives students a true sense of delivering 100-120 compressions per minute.
- Lights continuously indicate rate progression, confirming that compressions are within the 100-120 cpm range recommended by industry guidelines.
- The yellow light will flash if the upper limit of 120 cpm is exceeded, warning students to slow down.
- Instructors can monitor multiple students quickly and easily.
- In addition to the visual CPR rate monitor, students will also hear a clicking sound when the chest is pushed to the appropriate depth of 5 to 6 centimetres, allowing them to experience the actual force required to deliver chest compressions in real life.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- The CPR rate monitor can be purchased separately to upgrade any professional, unsupervised PRESTAN manikin.

INCLUDE:

- 2 Adult manikin Prestan
- 1 Child manikin Prestan
- 5 Baby manikin On loan
- 10 Lung replacements for adults Lend
- 5 Child Lung Replacements Loaned
- 10 Lung replacements Prestan Baby
- Instruction manual
- Carrying bag

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>