



STANDARD BREAST EXAMINATION TORSO

PRODUCT NUMBER: 1140



DESCRIPTION:

The **New Breast Examination Simulator** combines state-of-art materials to create a realistic look, feel and texture in addition to lifelike softness and durability. The breasts, nipples and areola have **realistic features**. The **cysts and masses located in the breasts** have distinct textures useful in recognizing normal and abnormal conditions found in breasts.

As such, the simulator provides an excellent platform on which trainees can **learn how to perform a breast examination** as well as demonstrate breast self-examination to women who want to learn how to do it.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



SKILLS:

- Learn to perform breast examination competently (i.e. spiral and up-and-down method of palpation)
- Breast awareness and self-examination

CASES AND PATHOLOGIES:

- Enlarged lymph node
- Fibroadenoma
- Fibrocystic mass
- Cyst
- Various sized masses

CHARACTERISTICS:

FEATURES:

- Upper torso with detachable breasts
- Left breast contains four masses in the breast tissue and two masses in the axilla region. The masses range in size from 14 to 19 mm, and their depths range from 6 to 16 mm beneath the surface.
- Right breast contains enlarged lymph node, fibroadenoma, fibrocystic mass and fluid-filled cyst (ranging from 12 to 24 mm)
- Breasts are attached to a manikin torso and can be easily removed and reassembled
- Simulator can be used in either the upright or reclining position
- Realistic texture and look

RIGHT BREAST LOCATION OF MASSES AND CYSTS:

- Enlarged lymph node – Diameter: 13mm Depth: 9mm
- Fibroadenoma – Diameter: 12mm Depth: 11mm
- Fibrocystic mass – Diameter: 24mm Depth: 8mm
- Cyst – Diameter: 19mm Depth: 6mm

LEFT BREAST LOCATION OF MASSES AND CYSTS:

- Mass 1 – Diameter: 16mm Depth: 13mm
- Mass 2 – Diameter: 15mm Depth: 6mm
- Mass 3 – Diameter: 18mm Depth: 6mm

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Mass 4 – Diameter: 17mm Depth: 16mm
- Mass 5 – Diameter: 14mm Depth: 10mm
- Mass 6 – Diameter: 19mm Depth: 13mm

INCLUDE:

- 1 Upper torso
- 1 Right breast containing 2 masses and 2 fluid filled cysts
- 1 Left breast containing 6 masses of varying sizes
- 1 CD with Breast Examination Learning Package
- 1 Soft carrying bag
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>