


# NEWBORN ARM FOR ULTRASOUND-GUIDED IV TRAINING

PRODUCT NUMBER: 2750



## DESCRIPTION:

**Baby Arm** is an anatomically correct, highly realistic arm phantom that can be used for ultrasound-guided IV training for nurses and physicians.

The phantom is an ideal investment to learn the skills of giving an intravenous injection on a baby's arm under ultrasound guidance.

## SKILLS:

- Ultrasound-guided venous access skills training
- Probe manipulation
- Recognition of anatomical structures

## CHARACTERISTICS:

### FEATURES:

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Mimics like a baby arm
- Great skin quality
- Anatomically correct Blood vessels
- Compatible with Ultrasound imaging
- Tested for IV in our labs
- The baby arm phantom is specifically designed for ultrasound-guided needle placement without involving the blood drawing procedure.
- The phantom is an ideal investment to learn the skills of giving an intravenous injection on a baby's arm under ultrasound guidance and it is made from durable and long-lasting materials.
- The blood vessels are prefilled with an ultra-soft rubber-like material so the phantom will never leak or require any extra blood mimicking fluid.

## INCLUDE:

- 1 Ultrasound-compatible newborn arm trainer
- 1 Box

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>