

# PEDI - Advanced pediatric care - 1 year

PRODUCT NUMBER: 1283



## DESCRIPTION:

The One-Year-Old PEDI® Nursing Care Patient Simulator allows participants to practice and develop neonatal nursing care skills. Includes an intubatable airway, injection training arm, and an intraosseous injection leg so that participants can work hands-on during chest compression and ventilation exercises

## SKILLS:

- Bandaging
- Hygiene care
- Palpation of pulses
- Ear canal examination
- BVM ventilation

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

- Airway management
- Nasal intubation
- Intubation, suctioning and Sellick maneuver
- NG and OG tube feeding
- Gastric suction
- Intramuscular injection exercises
- Intraosseous infusion
- IV training arm with intramuscular, intradermal and subcutaneous

## CHARACTERISTICS:

### FEATURES:

- Full-body with soft, lifelike face skin and molded hair
- Eyes can be closed or opened manually
- Head tilt/chin lift
- Endotracheal intubation
- Nasogastric intubation
- Tongue and teeth
- Intraosseous infusion
- IV training arm
- Femoral vein
- Manual palpable pulses: femoral and brachial
- BVM with realistic chest rise
- CPR anatomical landmarks
- Lavage/gavage
- Intramuscular injection sites
- Simulated ear canal
- Practice gastric suctioning exercises
- Articulating head, arms, and legs
- Realistic hands, feet, fingers, toes

## INCLUDE:

- 1 Full-size pediatric manikin
- 1 Set of clothes (T-shirt + shorts)
- 1 Soft carrying bag
- 1 Instructions manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>