



LOWER STUMP BANDAGING SIMULATOR

PRODUCT NUMBER: 1464



DESCRIPTION:

Practice and demonstrate **standard bandaging procedures**.

Can be used to demonstrate the attachment of prosthetic devices.

CHARACTERISTICS:

- Dimensions: 89 x 38 x 25 cm
- The lower torso is long enough to allow carrying the bandage around the body over the uninvolved hip at the level of the iliac crest.
- Note: There is no embedded, palpable skeletal material.
- Can be used to demonstrate the attachment of prosthetic devices.
- The lower torso is representative of a patient in a supine position with legs slightly abducted.
- One leg is amputated below the knee the other at mid-thigh.

INCLUDE:

- 1 Lower stump bandaging simulator
- Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>