


# AirSim Pierre Robin X

PRODUCT NUMBER: 1102



## DESCRIPTION:

TruCorp AirSim Pierre Robin X was developed in collaboration with clinicians from Children's Hospital Los Angeles. This realistic training manikin highlights difficult airway management techniques for infants with Pierre Robin Sequence/Syndrome (PRS). AirSim® X airway with 5-year warranty certified for over 20,000 intubation cycles.

## SKILLS:

- Oral and nasal intubation
- Bag/mask ventilation
- Laryngeal Mask Airway (LMA) and supraglottic device insertion
- Combitube insertion

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Double lumen endotracheal tube insertion
- Diagnostic bronchoscopy
- Lung isolation techniques using left and right endobronchial tubes and bronchial blockers
- Lung suctioning technique

## CHARACTERISTICS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- Illustrates the various congenital defects of an infant patient with Pierre Robin Sequence including; significant mandibular hypoplasia, glossoptosis, cleft palate and a bifid uvula
- User feedback through the visualization of stomach inflation and lung expansion
- Modular design allowing for all parts to be changed within your organization (if required)
- Delivered ready to use in a durable carry case for easy transportation and safe storage

## INCLUDE:

- 1 AirSim Pierre Robin X (assembled and ready to use)
- 1 TruCorp carry case
- 1 USB pen drive containing product user manual
- 1 bottle of TruCorp lubrication (100ml)

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>