




# SUPER OB SUSIE - Ultra-realistic childbirth simulator

PRODUCT NUMBER: 1141



## DESCRIPTION:

The new **SUPER OB SUSIE** is a **realistic childbirth torso** designed to aid educators **teaching labor and delivery management skills to learners of all levels**. The anatomically accurate structures allow learners to develop clinical skills transferable to scenario-based training. The modular design allows the functionality to expand with growing specific training needs.

To further extend use, additional modules of Post-partum hemorrhage and episiotomy repair are available.

## SKILLS:

- Execute cephalic deliveries

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Perform Ritgen's maneuver
- Practice breech deliveries (complete, frank, and footling)
- Perform Pinard's maneuver
- Simulate shoulder dystocia
- Practice emergency interventions: McRobert's maneuver, suprapubic pressure, posterior arm sweep and Wood's Screw
- Perform placental and umbilical exercises: placental delivery, partial placenta previa, nuchal cord, cord prolapse, true knot, and many others
- Practice Umbilical cord procedures like clamping and cutting
- Demonstrate internal rotation, expulsion, and external rotation
- Perform operative deliveries with real instruments like vacuum suction delivery and forceps assisted delivery
- Simulate nasal and oral suctioning
- Rectal suppository administration.

## CHARACTERISTICS:

### GENERAL:

- Adult-sized lower torso from diaphragm to quadriceps
- Smooth skin with seamless joints for enhanced realism
- Anatomically accurate pelvic landmarks including ischial spines and coccyx
- Supports bench top training or hybrid simulation
- Removable pregnant abdominal cover
- Articulating hip joints aid in performing McRobert's maneuver
- Removable cervix and birth canal
- Lifelike placenta with removable cotyledons
- Practice catheterization exercises with patent urethra
- Urine reservoir holds up to 400m
- Patent rectum accepts suppositories

### BIRTHING BABY:

- Full-term neonate has palpable landmarks, articulating limbs, and soft full-body skin
- Seamless, articulated joints
- Nasal and oral cavities
- Palpable fontanelle and sutures
- Smooth, seamless scalp
- Realistically articulating neck and torso
- Detachable umbilical cord and cord stumps for additional cutting and clamping

## INCLUDE:

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- 1 Super OB Susie Childbirth training torso
- 1 Pregnant abdomen
- 1 Full term neonate
- 1 Birthing handle
- 1 Placenta with cotyledons
- 2 Umbilical cord
- 2 Umbilical stump
- 2 Cervix
- 2 Birth canal
- 1 Talcum powder
- 1 Mineral oil
- 1 Urinary fill kit
- 1 Instructions manual

## ADDITIONAL MODULES:

### **POST-PARTUM HEMORRHAGE KIT (sold separately):**

- The PPH management package allows for postpartum skills training and emergency management.
- Includes: PPH perineum, blurred uterus, PPH cervix and blood concentrate.

### **EPISIOTOMY AND POST-PARTUM SUTURE KIT (sold separately):**

- Episiotomy repair trainers simulate human tissue and support the use of real sutures and instruments.
- Includes: Left Mediolateral Episiotomy with Perineum and Midline Episiotomy with Perineum.

### **BIRTH ASSIST MECHANISM (sold separately):**

- The manual lever of the delivery assist mechanism makes descent and rotation of the control effortless.
- Includes: Delivery mechanism, back panel and fixed fixings.

## **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

## **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>