

LEG TRAINER FOR HEMORRHAGE AND BLEEDING CONTROL

PRODUCT NUMBER: 2517



DESCRIPTION:

Our freestanding **trauma lower leg** is designed for first aid and military training and hemorrhage control.

This model can be customized for your specific needs.

SKILLS:

- Hemorrhage and bleeding control training
- Tourniquet application
- Usage of hemostatic agents

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Bandaging procedures training

CHARACTERISTICS:

GENERAL:

- Natural size leg trainer
- Integrated bleeding tubes
- Skin safe, prosthetic grade silicone
- Weight: 7 kg

FEATURES:

- This Full Leg Major Injury Task Trainer has an integrated major injury to the lower leg, exposing internal bone and muscle.
- Featuring integrated bleeding tubes, designed to provide realistic haemorrhagic bleeding using a manual pump, this task trainer enables students to practice haemorrhage control using tourniquets, bandages and haemostatic agents.

CUSTOMIZED WOUNDS (choose one of the following):

- Deep explosion wound with or without bone fragments
- Deep cut wound with embedded shrapnel
- Open Tib/Fib fracture
- Wide shot wound with lots of embedded pellets
- Deep cut wound from chainsaw
- Busted knee with a deep cut, grazes, and bruises

INCLUDE:

- 1 Freestanding silicone Leg model with embedded bloodlines
- 1 Customized wound scenario
- 1 Storage case
- Instructions manual

ADDITIONAL MODULES:

SIMULATED CLOTTING BLOOD (sold separately):

- The most realistic blood simulant currently available! Developed for training first responders in traumatic hemorrhage control, wound packing, and the use of advanced hemostatic agents.
- Use in conjunction with our Z fold or rolled Hemostatic Gauzes for training to create realistic clots.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Concentrate gel can be prepared and stored until needed or easily transported to the training area, where it can be reconstituted on site.
- The simulant accurately matches the color, opacity, and flow characteristics of human blood. It is nonbiological and nonhazardous, washes from skin and clothing with soap and water.

HEMOSTATIC Z-FOLDED GAUZE (sold separately):

- Use this simulated Hemostatic Gauze in conjunction with our simulated blood to initiate the simulated clotting.
- For training purposes only.
- Each package contains 3 individual bags. Each bag contains one strip of gauze measuring 7.5cm x 3.6m

HEMOSTATIC ROLLED GAUZE (sold separately):

- Use this simulated Hemostatic Gauze in conjunction with our simulated blood to initiate the simulated clotting.
- For training purposes only.
- Each package contains 3 individual bags. Each bag contains one strip of gauze measuring 7.5cm x 3.6m.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>