


# AirSim Child X

PRODUCT NUMBER: 1103



## DESCRIPTION:

The AirSim® Child X was developed in conjunction with the Children's Hospital Los Angeles. This anatomically correct manikin is based on true CT data from a 6 year old patient and includes the AirSim® X airway with 5-year warranty. The results are both realistic and anatomical pediatric features perfect for practicing basic airway management skills.

## SKILLS:

- AirSim®; Pediatric Airway with 5-year warranty
- Realistic feedback during bag/mask ventilation, supraglottic airway insertion and both naso- and orotracheal intubation
- Lifelike external and internal anatomical features allow for easy and successful demonstration of all pediatric video laryngoscopes (VL)

## MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

## CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Inflatable tongue bulb to simulate tongue edema
- Missing front teeth resemble natural dentals of a six year old child
- 'Real feel' skin covering for added realism during airway management procedure training

## CHARACTERISTICS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- Inflatable tongue bulb to create tongue edema providing a variance during training
- User feedback through the visualization of stomach inflation and lung expansion
- Missing front tooth to illustrate realistic dentals of a 6-year old
- Modular design allowing for all parts to be changed within your organization (if required)
- Can be upgraded to an AirSim Child Bronchi X for bronchoscopy training techniques
- Delivered ready to use in a durable carry case for easy transportation and safe storage

## INCLUDE:

- 1 AirSimChildX (assembled and ready to use)
- 1 TruCorp carry case
- 1 USB pen drive containing product user manual
- 1 bottle of TruCorp lubrication (100ml)

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>