




ULTRASOUND BREAST PHANTOM SIMULATOR

PRODUCT NUMBER: 1157



DESCRIPTION:

The **Ultrasound Breast Phantom Simulator** features a left and right breast attached to an adult upper torso to help participants practice and develop essential breast examination skills. The **left breast permits ultrasound identification of cysts versus dense masses**, while the right breast permits identification of cysts of different sizes and depths.

SKILLS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Learn to do ultrasound-guided needle aspiration
- Perform breast ultrasound imaging with real equipment
- Breast awareness and self-examination (note: the breast masses and cysts can only be located by ultrasound).

CASES AND PATHOLOGIES:

RIGHT BREAST LOCATION OF MASSES AND CYSTS:

- Cyst 1 – Diameter: 11 mm
- Cyst 2 – Diameter: 12 mm
- Cyst 3 – Diameter: 16 mm
- Cyst 4 – Diameter: 14 mm
- Cyst 5 – Diameter: 10 mm
- Cyst 6 – Diameter: 13 mm
- Cyst 7 – Diameter: 12 mm
- Cyst 8 – Diameter: 10 mm
- Cyst 9 – Diameter: 8 mm
- Cyst 10 – Diameter: 12 mm

RIGHT BREAST LOCATION OF MASSES AND CYSTS:

- Cyst 1 – Diameter: 14 mm
- Solid Mass 1 – Diameter: 10 mm
- Cyst 2 – Diameter: 10 mm
- Solid Mass 2 – Diameter: 15 mm
- Solid Mass 3 – Diameter: 20 mm
- Solid Mass 4 – Diameter: 13 mm
- Cyst 3 – Diameter: 8 mm
- Solid Mass 5 – Diameter: 16 mm
- Solid Mass 6 – Diameter: 15 mm

CHARACTERISTICS:

FEATURES:

- Visualize cysts and masses using real ultrasound equipment
- Learn to do ultrasound-guided needle aspiration in a relaxed environment before moving onto actual patients
- Left breast contains three cysts and six masses randomly placed in the breast and axilla
- Right breast contains ten cysts of different sizes and depths
- Remove fluid while seeing the cyst shrink and then disappear

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Realistic texture and look
- Self-healing skin that leaves no mark after removal of aspiration or biopsy needles
- Breasts easily removed and assembled
- Use in either the upright or reclining position
- Medium skin tone standard. Light or dark skin tone optional at no extra cost
- Indiscernible masses and cysts prompt students to resort to ultrasound imaging

INCLUDE:

- 1 Upper torso
- 1 Right breast containing 10 soft, fluid-filled cyst of varying sizes
- 1 Left breast containing 6 dense masses and 3 cysts randomly placed
- 1 Soft carrying bag
- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>