


Ultra-realistic 8 year old child simulator - LifeCast

PRODUCT NUMBER: 3604



DESCRIPTION:

Ultra-realistic 8 year old child simulator, with realistic details, such as veins and underlying structures, hair and an anatomically accurate airway for the correct training of procedures. Created to reflect the challenges posed by children of this age, and the procedures associated with emergencies.

Created with a manageable weight and average adult height. They do not require electrical power to operate and are therefore ideal for on-site simulation sessions; bringing a new level of realism to medical training, they help to create emotional attachment and commitment to clinical practice-based training.

SKILLS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Realistic airway management
- Oral and nasal intubation
- Cardiopulmonary resuscitation (CPR)
- Needle decompression in 2nd and 5th intercostal space
- Advanced Life Support (ALS)

CHARACTERISTICS:

GENERAL:

- Representation of an 8-year-old girl
- Measurements: 132 cm and 24kg
- Can be configured in a variety of ethnicities
- Created from real people through the processes of life-casting and 3D scanning
- Careful detailing to ensure immersion
- Mobility for handling
- No electrical connections required

INCLUDE:

Ultra-realistic 8 year old girl simulator
Instruction Manual
Transport bag

ADDITIONAL MODULES:

- Bilateral intraosseous infusion in tibia (Ref. 3608)
- Immersion capacity
- Patient Monitor (Ref. 1010)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>