

Adult nasogastric catheterisation and tracheostomy care simulator

PRODUCT NUMBER: 1465



DESCRIPTION:

The adult simulator for nasogastric catheterisation and tracheostomy care has been developed to assist in the instruction, training and practice of vital tracheostomy skills and the care of patients with respiratory conditions.

For exceptional realism, hydroxyethylcellulose can be mixed with water to simulate the mucous fluids of a real patient.

SKILLS:

- Gastrointestinal care procedures via nasal and oral access.
- Practice dressing changes and cuff inflation.
- Practice intubation with a laryngoscope blade.

CHARACTERISTICS:

GENERAL:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- This simulator features realistic landmarks, trachea, oesophagus, lungs and stomach.
- Other realistic procedures include:
- Nasogastric and feeding tube care.
- Insertion
- Irrigation
- Gastric lavage removal
- Catheter insertion
- Nasoenteric and oesophageal tube care
- Oropharyngeal and nasopharyngeal care
- Tracheostomy tube care
- Fluid can be added to the lungs and stomach for realistic tracheostomy and suctioning.
- Allows the use of size 6 tracheostomy tubes (LF01170U and XLF01170U) and size 6 nasogastric tube.
- Dimensions: 25x40x68 cm

INCLUDE:

- Adult simulator for nasogastric catheterisation and tracheostomy care
- Spray pump lubricant
- One hydroxyethylcellulose packet
- Instruction manual
- Hard carrying case.

*Simulator does not include endotracheal tracheostomy tube or nasogastric tube.

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>