

# AirSim Child Combo Bronchi X

PRODUCT NUMBER: 1108



## DESCRIPTION:

TruCorp AirSim Child Combo Bronchi X features the uniquely constructed AirSim® X airway with 5-year warranty enhanced with palpable cricoid landmarks, laryngeal cartilages and tracheal rings. This versatile pediatric intubation manikin is ideal for practicing percutaneous tracheostomy, ventilation, bronchoscopy and needle & surgical cricothyrotomy.

## SKILLS:

- Oral and nasal intubation
- Needle and surgical cricothyrotomy
- Percutaneous tracheostomy
- Bag/mask ventilation
- Laryngeal Mask Airway (LMA) and supraglottic device insertion

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Combitube insertion; Double lumen endotracheal tube insertion
- Diagnostic bronchoscopy
- Lung isolation techniques using left and right endobronchial tubes and bronchial blockers
- Lung suctioning techniques

## CHARACTERISTICS:

- AirSim X airway is certified to practice 20,000+ intubation cycles without fail and is supported by a 5-year airway guarantee
- Realistic larynx which features palpable cricoid landmarks, laryngeal cartilages and tracheal rings to provide positive user feedback and easy identification of the sternal notch
- Features bronchi which provides internal anatomically correct features down to the fourth generation
- Cost-effective consumables which are quick and easy to interchange, perfect for a busy teaching environment
- Inflatable tongue bulb to create tongue edema providing a variance during training
- User feedback through the visualization of stomach inflation and lung expansion
- Optional subcutaneous fat tissue with blood can be added to demonstrate a more difficult tracheostomy or cricothyroidotomy procedure
- Missing front tooth to illustrate realistic dentals of a 6-year old
- Modular design allowing for all parts to be changed within your organization (if required)
- Delivered ready to use in a durable carry case for easy transportation and safe storage

## INCLUDE:

- 1 AirSimChildCombo Bronchi X (assembled and ready to use)
- 1 TruCorp carry case
- 1 USB pen drive containing product user manual
- 1 bottle of TruCorp lubrication (100ml)
- 5 larynx inserts (1 supplied on manikin & 4 spare)
- 1 wraparound neck skin (supplied on manikin)

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>