



# Highly realistic full-term newborn simulator - Baby Nina

PRODUCT NUMBER: 2054



## DESCRIPTION:

Baby Nina is a full Soft silicone Sim baby with a super realistic newborn look.

## SKILLS:

- Birthing and delivery exercises
- Normal or emergency C-section training
- Vaginal birthing in many models of birthing manikins
- Intubation and ventilation
- Blood draw from umbilical vein

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- New parenthood training; bathing and dressing the baby

## CHARACTERISTICS:

- Being a birthing baby he is flexible enough to fold into a pregnant belly or uterus and durable enough to withstand c-sections or vaginal birthing.
- The baby has lungs and can be intubated and ventilated with visible chest rise.
- Her mouth and nose can be cleaned out
- Her umbilical stump has 3 bloodlines and blood can be taken from it. (Blood needs to be inserted before each use.)
- Each baby is unique in colouring and has eyelashes and curly black hair.
- Since the hair is very vulnerable it is wetted down with silicone to create the newborn wet look
- She can be washed, bathed, dried and powdered
- Length: 50cm. Head circumference: 37cm
- Weight: 3,2kg

## INCLUDE:

- 1 Baby Nina
- 1 Set of clothes
- 1 Care kit: Baby powder, baby shampoo and dummy

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>