

Suture jejunum simulator

PRODUCT NUMBER: 3463



DESCRIPTION:

This two-layer small intestine offers the opportunity to train the surgical closure of the intestinal wall and to perform different types of anastomosis.

This two-layer intestine simulator is especially soft and flexible, perfectly imitating the characteristics of the delicate jejunal tissue of dogs and cats. The visual and textural demarcation between the tunica mucosa and the tunica muscularis offers the opportunity to perform various suturing practices.

SKILLS:

- Incisions
- Sutures

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Anastomoses

CHARACTERISTICS:

GENERAL:

- Reusable several times
- Tear-resistant and retains its shape and texture
- Practical handling
- Length: ca. 24cm
- Diameter: ca. 1.5 cm
- Layers:
 - Tunica mucosa
 - Tunica muscularis

INCLUDE:

- 1 Suture jejunum simulator

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>