

# CANINE IV LEG

PRODUCT NUMBER: 1954



## DESCRIPTION:

The trainer represents the left foreleg of a medium-sized dog in the sternal recumbent position, and features two interchangeable sleeves – normal fur for routine blood draws, and a “shaved” area for i.v. catheter insertion.

## SKILLS:

- Bandaging and restraining/assistance practice
- Blood draw exercises
- IV catheter insertion

## INCLUDE:

- 1 IV canine leg trainer
- 1 Normal fur sleeve
- 11 Shaved area sleeve
- 1 Storage box
- 1 Built-in IV pole
- 1 Blood powder

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

- 1 Syringe
- 1 IV catheter
- 2 Fluid supply bags



**MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

**CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>