

Feline bladder catheterisation simulator - Male

PRODUCT NUMBER: 3758



DESCRIPTION:

This model is an anatomically correct male feline that allows practice in the correct placement of the urinary catheter.

An inner casing can be removed from the cat, containing the bladder, pelvic bone and penis with a soft foam surround. The bladder can be filled with water from the filling tube at the back of the holder, to obtain real information and practice draining the bladder. The foreskin can be rolled back to expose the penis for probe insertion.

SKILLS:

- Suture practice
- Bladder catheterisation
- Bladder drainage practice

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- The model can be attached to a desk or table and has non-slip rubber feet to prevent it from moving during the procedure.
- Anatomically correct:
 - Bladder
 - Pelvic bone
 - Penis
- A replaceable suture block is placed on the back of the cat, with a mesh reinforcement in the silicone, which allows the catheter to be sutured to the skin. Possibility to practice suturing the catheter in place.
- Practice correct placement of the urinary catheter.
- Silicone penis and retractable foreskin.
- Bladder filled with water to obtain real information.
- Water tube is fixed at the back with a plug.
- Tabletop design with non-slip rubber feet.

INCLUDE:

- 1 male feline bladder catheterisation simulator
- Water filling tube
- Instruction manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>