



Adult auscultation torso - SAM 4 PLUS

PRODUCT NUMBER: 4399



DESCRIPTION:

The torso simulator for auscultation training brings new unique features to enhance the auscultation learning experience. The SAM4 Plus simulator enhances the teaching and learning process, providing greater flexibility for teachers. It provides versatility and adaptability according to your auscultation needs.

SKILLS:

- Auscultation techniques

CASES AND PATHOLOGIES:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Library of normal and abnormal heart, respiratory, bowel and murmur sounds/noises in the highest sound quality, authentically recorded from live patients:
- 37 heart sounds
- 11 cardiopulmonary combinations
- 26 lung sounds
- 9 bowel sounds
- 6 murmurs
- 12 echocardiogram videos
- 3 dynamic sounds

CHARACTERISTICS:

- Simulator control via smartphone or tablet
- SAM4 application, available for iOS and Android
- Intuitive interface streamlines the learning process
- Heart rate control to quickly adapt to different scenarios
- Allows the use of real clinical stethoscopes
- Bluetooth connection for sound
- Ability to share the training experience with large audiences by connecting to earpieces
- Includes a 360-degree rotating base that optimises the demonstration of auscultation points and the exact position of the simulator
- ECG waveforms for all heart sounds
- Palpable carotid pulse
- Available in 7 languages: English, French, German, Polish, Portuguese, Russian, Spanish and French
- LED lights at all auscultation points
- Weight: 11 kg

INCLUDE:

- SAM4 Portable Male Torso Auscultation System

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

- Power cable
- User's manual



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>