

Emergency caesarean section and delivery simulator - C-CELIA

PRODUCT NUMBER: 4383



DESCRIPTION:

Extremely realistic obstetric simulator for emergency caesarean section deliveries, allowing you to practice various types of caesarean section deliveries with realistic incisions and amniotic fluid. With C-Celia, you can teach the entire operation from the first incision to the extraction of the foetus and delivery, and up to the repair of the uterus and the closure of the abdominal skin. The material the simulator is made of offers the trainee a realistic, hands-on training experience (no hard plastics); it is also very durable and easy to clean.

It is essential that professionals are trained and prepared for any eventuality in the delivery room.

SKILLS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- Complete transverse caesarean section (Pfannenstiel incision)
- Vertical emergency caesarean section (from umbilicus to symphysis pubis)
- Practice of caesarean sections with incisions
- Emergency foetal extraction and clinical protocol procedures
- Practice of suturing techniques and realistic soft tissue repair of caesarean section tissue
- Postpartum haemorrhage management

- Procedures, repair of complex lacerations and uterine artery injuries

- Supports PPH management techniques

- Bakri Balloon
- B-Lynch suture
- Hayman suture
- O'Leary suture

CHARACTERISTICS:

GENERAL:

- Caesarean section practice with incisions through the abdominal wall, muscle, fascia, linea alba, uterus and amniotic sac.
- Realistic anatomy allows for realistic simulations.
- Manually recreate complex foetal presentations in order to simulate foetal extraction situations to be performed by obstetricians and the medical team.
- Practice the medical procedures and clinical protocols necessary to perform caesarean sections and learn to work together and master the situation.
- Caesarean section through the abdominal wall and uterus with incisions that pass through:
 - Abdominal wall
 - Muscle
 - Fascia
 - Linea alba
 - Uterus
 - Amniotic sac

- Teamwork and human factors training for emergency caesarean section

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



INCLUDE:

- Emergency caesarean section simulator with pregnant abdomen
- 2 pieces for uterine incision
- 1 set of two newborn twins
- Repair kit: amniotic sac, adhesive to repair uterus and reuse 5 times (can be used up to 10 times)
- 5 abdominal accessories
- Lubricating gel

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>