

# SUPER OB SUSIE - Postpartum hemorrhage kit

PRODUCT NUMBER: 1670



## DESCRIPTION:

Postpartum hemorrhage package enables the training of postpartum emergency management skills.

## SKILLS:

- Practice catheterization with the urethra that connects to a realistic urinary bladder.
- Adjust uterine tone with the pump attachment
- Perform fundal massage.
- Practice Bakri balloon tamponade placement.
- Rectal suppository administration.

## CHARACTERISTICS:

### FEATURES:

Postpartum hemorrhage capability is available through the blood reservoir that holds up to 400mL of fluid.

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

## INCLUDE:

- 1 Postpartum hemorraghe perineum insert
- 1 Boggy uterus
- 1 Postpartum hemorraghe cervix
- 1 Elevating pillow
- 1 Blood concentrate
- 1 Fill kit
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>