

CERVICAL EXAM AND PAP SMEAR TEST TRAINER

PRODUCT NUMBER: 2280



DESCRIPTION:

The Pap test is the most frequently performed cancer screening test in the world. It is essential for women's healthcare practitioners to achieve competence in this vital but highly sensitive procedure.

CASES AND PATHOLOGIES:

- Normal
- Early pregnancy
- Polyp
- Early cancer
- Late cancer
- Inflammation
- Dysplasia

CHARACTERISTICS:

- With the Cervical Exam and Pap Test Trainer, students can practice external and internal examinations, speculum insertion, and cytological

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



sampling techniques on a simulator that is incredibly realistic and easy to use.

- Ideal for patient education and simulation scenarios.
- Disassembles completely for easy cleaning
- Extra-soft material and lifelike detail give extremely realistic training experience
- Fraction of the cost of a full-size simulator
- Includes materials for making realistic body fluids and secretions
- Modular design offers quick change-out for different scenarios
- Movable cervix
- Compact size for easy setup, storage, and transport.

INCLUDE:

- External genital pad with vagina
- 7 Cervices
- Support blocks
- Baby powder
- Hydroxyethyl cellulose
- Blood packet
- Syringe
- Lubricant
- Instruction manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>