

# MIGS SC - Suprachoroidal migs

PRODUCT NUMBER: 1814



## DESCRIPTION:

The **MIGS-SC model** features realistic angle visualization and enables simulated **implantation of piercing type MIGS** (istent, xen) and **supra-choroidal MIGS** (Supra, CyPass, miniject).

Moreover, you can develop **surgical gonioscopy skills** and practice the bi-manual technique.

## SKILLS:

- Surgical gonioscopy skills
- Bi-manual technique
- Angle visualization
- Insertion of MIGS devices into the supra-choroidal space

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Instrument manoeuvring skills

## CHARACTERISTICS:

These eyes have been designed for the insertion of a wide variety of MIGS devices and to practice gonioscopy skills.

### INSTRUCTIONS FOR CARE:

Follow these recommendations to maximize the life of your models:

- Store in a cool, dry and darkplace (a drawer will be fine). Extended exposure to some indoor lights or sunlight (UV) may affect material properties. Prolonged exposure to humidity or high temperatures may adversely affect material properties.
- Do not place heavy objects on top of the model's box. Prolonged compression may deform the models.

## INCLUDE:

- 1 Supra-choroidal MIGS eye simulator
- 1 Storage box

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>