

# NEWBORN IV AND ARTERIAL ACCESS ARM

PRODUCT NUMBER: 1933



## DESCRIPTION:

The new Newborn IV And Arterial Access Training Arm is a realistic and versatile task trainer designed for teaching intravenous and arterial access skills to learners of all levels.

## SKILLS:

- Palpation of radial and brachial pulse
- Venous access exercises
- Venous and arterial cannulation
- Infusion exercises
- Blood draw and collection
- Arterial blood sampling and collection

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



## CHARACTERISTICS:

### GENERAL:

- Lifelike skin provides a realistic look and feel
- Realistic arterial and venous network
- Adjustable pulse strength and rate (0 - 180 BPM)
- External blood reservoir for continuous fluid supply and vessel pressurization and observe flashback
- Easily replaceable veins and arteries
- Articulating wrist allowing palmar flexion
- Use as a stand-alone or attached to stand
- Allows transillumination for vein access
- Latex-Free

### ARTERIAL AND VENOUS INJECTION:

- Ulnar and cubital venous access sites support cannulation, infusion, bolus injection, blood draw/collection
- Radial and brachial arterial access sites support cannulation and blood sampling/collection
- Palpable radial and brachial pulse

## INCLUDE:

- 1 Newborn IV And Arterial Access Training Arm
- 1 Stand
- 1 Blood bag reservoir
- 1 Control box
- 1 Blood dispensing bag
- 1 Accessories set
- 1 Power supply
- 1 User manual

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>