

Arterial blood gas simulator

PRODUCT NUMBER: 4371



DESCRIPTION:

This arterial blood gas simulator is ideal for practice and demonstration in arterial blood sampling for blood gas analysis. It allows confirmation of correct needle placement in the artery.

SKILLS:

- Identification of suitable sites for arterial puncture
- Palpation of the pulse
- Arterial blood sampling
- Blood gas analysis

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- Realism of appearance and execution necessary to provide realistic practice for trainees.
- Automatic pressure produces a backflow of blood into the syringe after correct puncture.
- Hand skin and artery can be replaced when necessary to change them
- The simulator is latex-free
- Dimensions: 60 x 23.5 x 18.5 cm approx.
- Weight: 5.5 kg

INCLUDE:

- Arterial blood gas simulator

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>