

Adult torso for intraosseous access and sternal infusion training

PRODUCT NUMBER: 2737



DESCRIPTION:

The Adult Intraosseous Access and Sternal Infusion Simulator is the first of its kind for dramatic training designed to demonstrate, simulate and prepare all levels of healthcare professionals for the use of the intraosseous sternal approach.

Created with the guidance of experts who developed this technique, the incredible realism and accurate landmarks will allow trainees to simulate this new life-saving procedure with skill, confidence and precision. This procedure is a quick and easy alternative to conventional intravenous infusion. It saves time and saves lives in demanding emergency settings.

SKILLS:

- Intraosseous access
- Sternal infusion

CHARACTERISTICS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- The simulator uses a replaceable bone designed to be rotated and punctured several times.
- The sternum can be filled with simulated blood to better visualise the actual flashback.

INCLUDE:

- Training simulator
- 10 replaceable simulated bones
- Syringe with tube
- Blood mixture
- 1 hard carrying case.