

# PEDIATRIC INJECTION LEG - 1 YEAR

PRODUCT NUMBER: 1169



## DESCRIPTION:

The Mike® and Michelle® Intraosseous leg simulates that of an average 1 year-old child. The leg is attached to a half section of a lower torso. A replaceable tibial bone lies under a smooth outer skin, and is molded with anatomic landmarks for teaching intraosseous access and infusion. A femoral vein/ artery pair and an intramuscular injection site are included.

## SKILLS:

- Intraosseous access
- Infusion of medicaments/fluid exercises
- Blood draw exercises
- Intramuscular injection
- IV training exercises

## CHARACTERISTICS:

### FEATURES:

- Realistic "pop" as needle penetrates cortex and enters bone marrow

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- Sixteen (16) interchangeable, realistic tibial bones
- Blood return upon successful entry
- Infusion of fluid/medications
- Easily replaceable skins and veins that withstand multiple needle sticks
- Mounted on smoked Lucite® base and supplied with blood bag and stand, synthetic blood concentrate, squeeze bulb, and spare leg skin

## INCLUDE:

- 1 IO leg with stand
- 16 Interchangeable tibial bones
- 1 Blood bag
- 1 Squeeze bulb
- 1 Spare leg skin
- 1 Funnel
- 1 Blood concentrate
- 1 Soft carrying bag
- 1 Instructions manual

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>