

PETER PICC LINE - UPPER TORSO FOR PICC

PRODUCT NUMBER: 2093



DESCRIPTION:

Peter PICC Line™ permits the practice of placing Peripherally Inserted Central Catheter Lines (PICC Lines).

The rapid acceptance and increased use of PICC lines creates a need to be able to properly teach this procedure. Peter PICC Line™ allows the teaching, practice and assessment of this insertion technique.

SKILLS:

- PICC insertion exercises and assessment
- Needle placement
- Measurement and correct placement of catheters
- Sterile technique skills training

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- Visualization of underlying cephalic, basilic and median basilic veins

CHARACTERISTICS:

- Peter PICC Line™ permits the practice of placing Peripherally Inserted Central Catheter Lines (PICC Lines).
- The rapid acceptance and increased use of PICC lines creates a need to be able to properly teach this procedure.
- Peter PICC Line™ allows the teaching, practice and assessment of this insertion technique.
- This portable teaching model has a replaceable translucent arm skin that allows the visualization of the underlying cephalic, basilic and median basilic veins which are presented as barely visible blue veins.
- The superior vena cava, subclavian, axillary and jugular veins are also present, plus a movable chin that simulates the occlusion of the jugular to prevent the PICC Line from traveling this route.
- Proper placement can be confirmed by viewing the distal tip of the PICC in the viewable superior vena cava on the side of the model.
- Comes with removable ribs and overlying tissue pad.
- The ribs are palpable permitting the practice of measuring catheter length to the correct intercostal space and confirming proper measurement by viewing the distal tip in the SVC.

INCLUDE:

- 1 Body
- 2 Arm Vein
- 1 Body Vein
- 2 Arm skin
- 1 Body Skin
- 1 Ribs
- 1 Muscle Tissue
- 1 SVC Viewing Vial
- 1 Large Reservoir bag
- 1 Pint simulated blood
- 1 Carrying Case

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>