


X-RAY PHANTOM PELVIS - TRANSPARENT

PRODUCT NUMBER: 1938



DESCRIPTION:

This X-Ray part phantom gives the unique opportunity to take x-ray images of a pelvis with 2 lumbar vertebrae and femur stumps again and again.

The Phantom includes real human bones and allows taking real x-ray images.

SKILLS:

- Positioning and alignment techniques in projection radiography
- Research applications

CHARACTERISTICS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- The Phantom includes real human bones and allows taking real x-ray images of a pelvis with 2 lumbar vertebrae and femur stumps.
- The model is perfect for schools and education, but also for medical technicians since the same bones can be x-rayed again and again in different settings without the danger of harming a patient.
- The bones are embedded in transparent plastic.
- This phantom is handmade and unique. It may differ in size and shape.

INCLUDE:

- 1 Pelvis phantom

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>