

Episiotomy and perineal suture trainer

PRODUCT NUMBER: 1153



DESCRIPTION:

Simulated birth assistance training to improve clinical skills: The episiotomy and suture trainer allows for practice in the professional care and suturing of vaginal birth wounds.

It is recommended that medical personnel (midwives and doctors) receive continuous training in the diagnosis and treatment of birth injuries.

In order to realistically reproduce these birth wounds, the simulator allows the following muscle structures to be palpated and sutured in the perineal area:

- M. sphincter ani externus (external sphincter ani externus muscle)
- M. transversus perinei superficialis (superficial transverse muscle of the perineum)
- M. bulbospongiosus (bulbospongiosus muscle)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

SKILLS:

This trainer enables wound management to be demonstrated:

- Median episiotomy
- Mediolateral episiotomy
- Lateral episiotomy
- Vaginal lacerations
- Labial lacerations
- Perineal suture

CHARACTERISTICS:

GENERAL:

- Muscle and skin layers can be sutured separately.
- Material resistant to multiple cuts, repetition.
- Realistic fabrics to the touch.
- Easy care and cleaning.
- Compact trainer.
- The adjustable support of the simulated tissues is stable.
- Detachable parts, you can work without the support.
- Silicone-based and latex-free material.

INCLUDE:

- Episiotomy and Perineal Suture Trainer – Silicone Part
- Base
- Support

ADDITIONAL MODULES:

OPTIONAL MODULE:

- Refill of episiotomy pads (5pcs) – Ref. 3713

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>