



Colon simulator with suture mesentery

PRODUCT NUMBER: 3377



DESCRIPTION:

In addition to a mucous membrane and a muscularis tunica, this intestinal simulator with attached mesentery offers a particularly realistic representation of the various intestinal layers. In this way, incisions, intestinal sutures and various anastomoses can be performed under realistic conditions. The tunica mucosa, tunica muscularis and mesentery are optically separated from each other, so that the application of various sutures can be trained.

SKILLS:

- Incisions
- Intestinal sutures
- Anastomosis

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- Reusable several times
- Tear-resistant and retain their shape and texture
- Practical handling
- Length: ca. 24 cm
- Diameter: ca. 2.7 cm
- Layers:
 - Mucous tunic
 - Muscular tunic
 - Mesentery

INCLUDE:

- 1 Colon simulator with suture mesentery

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
Spain

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>