

# DIAGNOSTIC AND PROCEDURAL EAR TRAINER

PRODUCT NUMBER: 2155



## DESCRIPTION:

Now with this new trainer, students, residents, and practicing physicians can **master the skills needed to examine the human ear** using visual cues, correctly **diagnose common diseases, clean the ear canal, remove a foreign body**, and perform a **myringotomy** with ear tube insertion.

## SKILLS:

- Examination of the ear
- Diagnosis of normal and abnormal conditions
- Earwax removal exercises
- Removal of foreign bodies in the ear canal
- Myringotomy

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



## CASES AND PATHOLOGIES:

- Normal tympanic membrane
- Mucoid otitis media
- Serous otitis media with an air-fluid level
- Chronic otitis media with small and large tympanic membrane perforations
- Two views of attic cholesteatoma
- Atelectatic middle ear otitis
- Tympanic sclerosis

## CHARACTERISTICS:

- The Diagnostic and Procedural Ear Trainer is anatomically correct.
- Life-like right and left ears teach diagnostic and procedural techniques and are removable for easy maintenance and storage.
- For the ultimate in realism, the middle ear can be filled with various fluids.
- The instructor controls the color and consistency of the fluids.
- Another feature that brings this trainer as close to life as possible is the inclusion of nine diagnostic cartridges with full-colour photo prints
- Weight: 5kg

## INCLUDE:

- 1 Head model with two ears
- 1 Synthetic wax
- 1 Standard middle ear cartridge with syringe
- 1 Adjustable stand
- 1 Hard carry case

## ADDITIONAL MODULES:

- **PNEUMATIC OTOSCOPY KIT (option):**
- This option offers an additional feature for the Diagnostic and Procedural Ear Trainer [Ref.2155]
- While pneumatic otoscopy is indispensable for examination of the ear, studies have shown that most are performed incorrectly!
- Now with the valuable feedback offered by this kit, students, residents, and practicing physicians can learn how to generate the correct amount of air pressure needed to perform accurate pneumatic otoscopy to test for tympanic membrane mobility.
- Kit includes:

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



**More Than Simulators™**  
*Excellence. Every patient. Every time.*

- 1 Pneumatic pressure gauge
- 2 Pneumatic ears with pressure tubing
- 100 Eardrums membranes



**MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
Spain

**CONTACTO**

[sales@morethansimulators.com](mailto:sales@morethansimulators.com)  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>